

FILLOGY
FILLING EVOLUTION

FILLING KIT FOR
E-CIGARETTE LIQUIDS





FILLOGY 2+/250 + Filling set LIQUID M, Art. 301063 + 304185

Liquid and hot products, e.g.

- liquid cleaners
- food supplements for animals

Product feeding: via hose Ø 14 mm
 Filling speed: 7 l/min
 Filling amount: 50 ml - 1,000 ml
 Filling accuracy: +/- 1-2 %
 Filling temperature: 2 - 95 °C / 35,6 - 203 °F



FILLOGY LIQUID S AUTO S, Art. 301067

Liquid and flowable products, e.g.

- wood glues
- fine mechanics oils
- e-cigarette liquids

Product feeding: via hose Ø 14 mm
 Filling speed: up to 9 l/min
 Filling amount: 10 ml - 1,000 ml
 Filling accuracy: +/- 1-2 %
 Filling temperature: 2 - 95 °C / 35,6 - 203 °F



FILLOGY FLOWABLE M AUTO S, Art. 301071

Small amounts of flowable products, e.g.

- balms
- shoe polishes
- cleanser gels

Product feeding: via funnel 5 l
 Filling speed: up to 5 l/min
 Filling amount: 2 ml - 100 ml
 Filling accuracy: +/- 1-2 %
 Filling temperature: 2 - 95 °C / 35,6 - 203 °F



FILLOGY LIQUID S AUTO M, Art. 301024

Liquids products, e.g.

- glycerin
- molasse

Product feeding: via hose Ø 14 mm
 Filling speed: up to 9 l/min
 Filling amount: 500 ml - 2,000 ml
 Filling accuracy: +/- 1-2 %
 Filling temperature: 2 - 95 °C / 35,6 - 203 °F

